# D Zygosity-TYPE

BAGene product line

 **BAG** HEALTH CARE

IVD

CE

LOT 605ZY1

REF 6648

 10/2015

 2018-10

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*In presence of weak D type 4.0/4.1, 4.2, 29 and partial D type RHD-CE(8-9)-D false negative reactions in PCR mix no. 1 may occur.  
In case of RHD psi homo- or heterozygous false positive reactions in PCR mix no. 2 may occur.*

**Proben-ID / Sample ID:**

**Probenname / Name:**

**Geb.-Datum / Birthdate:**

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**Datum / Date:**

**Unterschrift / Signature:**

**Gelbild / Gel picture:**

Version 7.0 – 10/2016

BAG Health Care GmbH

Amtsgerichtsstraße 1-5  
35423 Lich/Germany

Tel.: +49 (0) 64 04/925 - 0  
Fax: +49 (0) 64 04/925 - 250

www.bag-healthcare.com  
info@bag-healthcare.com

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Tel.: +49 (0) 64 04/925 - 450  
Fax: +49 (0) 64 04/925 - 460  
verkauf@bag-healthcare.com

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## Short instructions

1. Take out the required number of plates/strips from -20...-80°C and thaw the 10 x PCR-buffer.
2. The first reaction mix can be identified by the printed lot number.
3. Pipet the mastermix, consisting of 10 x PCR-buffer, DNA-solution, Taq-Polymerase and distilled water and mix well. The different DNA-SSP Kits all work with the same mastermix and can therefore be combined.
4. After vortexing add **10 µl** of this mixture immediately to each dried reaction mixes.
5. Close the tubes tightly with the respective strip caps. Slightly move the plates/strips to dissolve the pellet at the bottom of the tube. The complete PCR-solution should settle on the bottom of the tube.
6. After PCR and separation of the amplicons in the gel, the evaluation can be performed.
7. Store all reagents at -20...-80°C after use.

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